

EZ-Dry 0.25 L Moisture Trap

For maintaining anhydrous conditions of solvents up to 0.25 L.

Contains: 1.0 g activated molecular sieves in an inert highly permeable porous membrane.

<p>Product Code: HR-0111</p>	
<p>Moisture Trap with molecular sieves</p>	
Quality Parameter	Specification
<p>Mass molecular sieves:</p> <p>Equilibrium water adsorption capacity, 50% relative humidity, 20°C, 24h:</p> <p>Residual Water Content, 550°C, 2h (molecular sieves only):</p>	<p>0.8 – 1.2 g</p> <p>≥ 22.0 %</p> <p>≤ 1.5 %</p>
<p>Notice:</p> <p>Storage:</p> <p>Shipping:</p>	<p>The use of this product is strictly limited to trained personnel for professional manufacturing, laboratory, or research purposes. Final Fitness-For-Use must be determined by and is the sole responsibility of the end-user.</p> <p>Material Safety Data Sheets are available upon request at info@empbiotech.com</p> <p>This product has been manufactured, tested and assembled under the control of an ISO 9001 registered quality system.</p> <p>Ambient.</p> <p>Ambient.</p>