

## EZ-Dry PRIME Moisture Trap, 20 L

For rapid reduction of water content in solvents up to 20 L.  
 Contains: 100 g activated Zeolite molecular sieves sealed in an inert highly permeable porous membrane.

<b>Product Code:</b> HR-0120	
Moisture Trap with Zeolite molecular sieves	
Quality Parameter	Specification
Mass Zeolite:	95 – 105 g
Equilibrium water adsorption capacity, 50% relative humidity, 20°C, 24h:	≥ 22.0 %
Residual Water Content, 550°C, 2h (zeolite only):	≤ 1.5 %
<p style="text-align: right;">Notice:</p> <p>The use of this product is strictly limited to trained personnel for professional manufacturing, laboratory, or research purposes. Final Fitness-For-Use must be determined by and is the sole responsibility of the end-user.</p> <p>Material Safety Data Sheets are available upon request at <a href="mailto:info@empbiotech.com">info@empbiotech.com</a></p> <p>This product has been manufactured, tested and assembled under the control of an ISO 9001 registered quality system.</p> <p style="text-align: right;">Storage:</p> <p style="text-align: right;">Shipping:</p>	<p>Ambient.</p> <p>Ambient.</p>