

EZ Dry ULTRA 50 g Moisture Trap

For maintaining anhydrous conditions of solvents up to 10 L

Contains: 50 g activated molecular sieves in an inert highly permeable porous membrane

<p>Product Code: HR-0137</p>	
<p>Moisture Trap with molecular sieves</p>	
Quality Parameter	Specification
<p>Mass molecular sieves:</p> <p>Equilibrium water adsorption test of molecular sieves</p> <p>Solvent Standard: 1600ppm water content, 20 °C</p>	<p>49 g – 51 g</p> <p>≤ 20 ppm water content after 24 h</p>
<p>Notice:</p> <p>Storage:</p> <p>Shipping:</p>	<p>The use of this product is strictly limited to trained personnel for professional manufacturing, laboratory, or research purposes. Final Fitness-For-Use must be determined by and is the sole responsibility of the end-user.</p> <p>Material Safety Data Sheets are available upon request at info@empbiotech.com</p> <p>This product has been manufactured, tested and assembled under the control of an ISO 9001 registered quality system.</p> <p>Ambient.</p> <p>Ambient.</p> <p>This product has passed all emp Biotech GmbH quality control specifications.</p>