

EZ Dry ULTRA 100 g Moisture Trap

For maintaining anhydrous conditions of solvents up to 20 L

Contains: 100 g activated molecular sieves in an inert highly permeable porous membrane

<p>Product Code: HR-0138</p>	
<p>Moisture Trap with molecular sieves</p>	
<p>Quality Parameter</p>	<p>Specification</p>
<p>Mass molecular sieves: Equilibrium water adsorption test of molecular sieves Solvent Standard: 1600 ppm water content, 20 °C</p>	<p>98 – 102 g ≤ 20 ppm water content after 24 h</p>
<p>Notice: Storage: Shipping:</p>	<p>The use of this product is strictly limited to trained personnel for professional manufacturing, laboratory, or research purposes. Final Fitness-For-Use must be determined by and is the sole responsibility of the end-user. Material Safety Data Sheets are available upon request at info@empbiotech.com This product has been manufactured, tested and assembled under the control of an ISO 9001 registered quality system.</p> <p>Ambient. Ambient.</p> <p>This product has passed all emp Biotech GmbH quality control specifications.</p>