

SMART Chromatography™ SP Boost SMART Chromatography™ IEX FPLC column with 10 mL bed volume and 12 cm bed height	
Product Code:	SC-4202-M010.0-012
Quality Parameter	Specification
Packing evaluation asymmetry: (no double peak, no shoulder)	$0.7 \leq A_s \leq 1.6$
Asymmetry LOT homogeneity:	$CV \leq 8 \%$
Notice:	<p>The use of this product is strictly limited to trained personnel for professional manufacturing, laboratory, or research purposes. Final Fitness-For-Use must be determined by and is the sole responsibility of the end-user.</p> <p>SMART chromatography™ is covered by US Patent 11,325,104 & 11,731,107.</p>